



The *smarter source* for affordable specialty instruments

MASTEL BAKEWELL BUBBLE LEVEL

PRE-OP MARKERS




- Toric/LRI marker, 8.2mm - 11.6mm marking fins, includes 3 extra levels
- Simplifies pre-marking a true horizon
- Creates peripheral marks at 0/180, and can be used with or without ink
- Simply hold the instrument like a dart, align the center of the marker with the patient's visual axis and align the marker fins to horizontal using the bubble
- Designed for Toric IOL's, LRI, AK, CRI and PCRI markings

**2-Ray
MRMT-BB**



Ambler Surgical is an authorized distributor for the Mastel Precision line of ophthalmic instruments and diamond scalpels. Mastel Precision has been involved exclusively in the design and production of ophthalmic instruments since 1981. Mastel Precision offers unparalleled design, quality and precision-crafted surgical instruments and diamond scalpels for astigmatism management and corneal procedures that are considered unmatched in the ophthalmic industry.

To view pricing and our complete line of Mastel Precision Instruments, visit us online.



OPHTHALMOLOGY

amblersurgical.com/mastel

1-888-407-0006