



The *smarter source* for affordable specialty instruments

MASTEL CIONNI DEGREE GAUGE RINGS

FIXATION/ALIGNMENT

Bores Handle



Thornton Handle




- 12.75mm Internal diameter / 16.00mm Outside diameter
- Ultra-Low Profile with Reverse Bevel
- Classic Mendez enumerations allow orientation to the steep corneal meridian
- Patterns (dots, arrows and diamonds) make it easier to determine arc length and orientation
- Designed for Toric IOL's, LRI, AK, CRI and PCRI markings

Handle Location	0/180°	90°
Bores Handle	RFBG-CR1216C-0	RFBG-CR1216C-9
Thornton Handle	RFTG-CR1216C-0	RFTG-CR1216C-9



Ambler Surgical is an authorized distributor for the Mastel Precision line of ophthalmic instruments and diamond scalpels. Mastel Precision has been involved exclusively in the design and production of ophthalmic instruments since 1981. Mastel Precision offers unparalleled design, quality and precision-crafted surgical instruments and diamond scalpels for astigmatism management and corneal procedures that are considered unmatched in the ophthalmic industry.

To view pricing and our complete line of Mastel Precision Instruments, visit us online.



OPHTHALMOLOGY

amblersurgical.com/mastel

1-888-407-0006