



The *smarter source* for affordable specialty instruments

## MASTEL GIMBEL-MENDEZ DEGREE GAUGE RINGS

### FIXATION/ALIGNMENT

Bores Handle



Thornton Handle




- 11.75mm Internal diameter / 16.00mm Outside diameter
- Bevel lends itself to use as a guide ring for LRI/AK when desired
- Ideal for corneal epithelial pre-marking approaches and has a higher profile in the periphery
- Designed for Toric IOL's, LRI, AK, CRI and PCRI markings
- Jarvi enumerated pattern available

| Handle Location | 0/180°      | 90°         | 135°        |
|-----------------|-------------|-------------|-------------|
| Bores Handle    | RFBG-GM1116 | RFBG-MG1116 | RFBG-KM1116 |
| Thornton Handle | RFTG-GM1116 | RFTG-MG1116 | RFTG-KM1116 |



Ambler Surgical is an authorized distributor for the Mastel Precision line of ophthalmic instruments and diamond scalpels. Mastel Precision has been involved exclusively in the design and production of ophthalmic instruments since 1981. Mastel Precision offers unparalleled design, quality and precision-crafted surgical instruments and diamond scalpels for astigmatism management and corneal procedures that are considered unmatched in the ophthalmic industry.

To view pricing and our complete line of Mastel Precision Instruments, visit us online.



OPHTHALMOLOGY

[amblersurgical.com/mastel](http://amblersurgical.com/mastel)

1-888-407-0006