




The *smarter source* for affordable specialty instruments

MICROSURGICAL ALUMINUM STERILIZATION TRAYS

- Aluminum, single-piece molded construction
- Highly durable, non-rusting and maintenance-free
- Dimpled base tray to avoid moisture under the mat
- All trays include base, lid, and silicone finger mat
(plus insert on double-level trays)

► SINGLE-LEVEL TRAYS



93-778 2 1/2" x 6 1/2" x 3/4", 2-4 Instruments

93-780 2 1/2" x 7 1/2" x 3/4", 2-4 Instruments

93-782 4" x 7 1/2" x 3/4", 6-8 Instruments

93-783 6 1/2" x 8 1/4" x 3/4", 10-12 Instruments

93-784 6" x 10" x 3/4", 12-15 Instruments


93-785 6" x 10" x 1 1/2" (deep), 12-15 Instruments

93-793 9" x 12" x 3/4", 18-22 Instruments

93-794 10" x 15" x 3/4", 25-35 Instruments

93-795 10" x 15" x 1 1/2" (deep), 25-35 Instruments

► DOUBLE-LEVEL TRAYS



93-786 6" x 10" x 1 1/2", 25-30 Instruments

93-787 6" x 10" x 3 1/4" (deep), 25-30 Instruments

93-796 10" x 15" x 1 1/2", 50-70 Instruments

To view pricing and our complete line of sterilization trays, visit us online.

