



The *smarter source* for affordable specialty instruments

## MASTEL DAVIS PLUMB BOB


### PRE-OP MARKERS

- Solomon plumb-bob 3-ray pre-op corneal gravity marker, 4 7/8", marks 3, 6, and 9 o'clock cardinal positions, Thornton titanium handle
- Davis plumb-bob 4-ray pre-op corneal gravity marker, 4 7/8", marks 3, 6, 9 and 12 o'clock cardinal positions, Thornton titanium handle
- 5.0mm Inside Diameter / 12.0mm Outside Diameter on horizontal
- 5.0mm Inside Diameter / 9.5mm Outside Diameter on vertical
- Markers pair well with the Zero Parallax Visual Axis Glasses for precise alignment on the visual axis



Ambler Surgical is an authorized distributor for the Mastel Precision line of ophthalmic instruments and diamond scalpels. Mastel Precision has been involved exclusively in the design and production of ophthalmic instruments since 1981. Mastel Precision offers unparalleled design, quality and precision-crafted surgical instruments and diamond scalpels for astigmatism management and corneal procedures that are considered unmatched in the ophthalmic industry.

To view pricing and our complete line of Mastel Precision Instruments, visit us online.



OPHTHALMOLOGY

[amblersurgical.com/mastel](http://amblersurgical.com/mastel)

1-888-407-0006