

Sterile Processing Department Essentials Catalog




Enhancing outcomes for
patients and their caregivers:

Hill-Rom

The following brands are registered trademarks of Aspen Surgical Products, Inc:
ColorTrack®, ErgoSupport®, Oxboro®, SafeClean®, SurgiSafe®, and Sterion®.

DuoGuard™ and InnerVent™ are trademarks of Aspen Surgical Products, Inc.

© 2016 Aspen Surgical Products, Inc. All rights reserved.




Oxboro® Instrument Tip Protectors

- Designed to protect valuable instruments from breaking and chipping during sterilization and handling
- Help eliminate reprocessing costs caused by torn pouches or wrap and protect staff from sharp sticks during handling

Osteotome Tip Protectors

- Help Protect Osteotome instrument edges from chipping with these tip protectors
- Available tinted or clear and vented or non-vented to comply with user preferences
- All types allow visual instrument inspection and are gas and steam permeable
- Not made with natural rubber latex and intended for single-use



Clear				Blue			
Part No.	Part No.	Size	Qty/Bag	Part No.	Part No.	Size	Qty/Bag
Vented	Non-Vented	mm		Vented	Non-Vented	mm	
094031BBG	094021BBG	6.4 x 19	50	094011BBG	094001BBG	6.4 x 19	50
094032BBG	094022BBG	9.5 x 25	50	094012BBG	094002BBG	9.5 x 25	50
094033BBG	094023BBG	12.7 x 25	50	094013BBG	094003BBG	12.7 x 25	50
094034BBG	094024BBG	19 x 25.4	50	094014BBG	094004BBG	19 x 25.4	50
094035BBG	094025BBG	25 x 25	50	094015BBG	094005BBG	25 x 25	50
094036BBG	094026BBG	32 x 32	50	094016BBG	094006BBG	32 x 32	50
094037BBG	094027BBG	38 x 38	50	094017BBG	094007BBG	38 x 38	50
094038BBG	094028BBG	44.5 x 32	50	094018BBG	094008BBG	44.5 x 32	50
094040BBG		Assorted	50	094020BBG		Assorted	50

Standard Tip Protector Sizing Guide

Use this guide to determine the standard tip protector on the following page that will best fit your instrument.



Standard Instrument Tip Protectors

- Available solid, tinted or clear and vented or non-vented
- Gas or steam permeable and can be used in flash or conventional autoclave processes
- Solid colors are radiopaque
- Not made with natural rubber latex and intended for single-use



1.6mm x 19mm	Part No.	Qty/Bag
Solid Vented	091011BBG	100
Solid Non-Vented	091001BBG	100
Tinted Vented	093011BBG	100
Tinted Non-Vented	093001BBG	100
Clear Vented	092011BBG	100
Clear Non-Vented	092001BBG	100



2mm x 19mm	Part No.	Qty/Bag
Solid Vented	091012BBG	100
Solid Non-Vented	091002BBG	100
Tinted Vented	093012BBG	100
Tinted Non-Vented	093002BBG	100
Clear Vented	092012BBG	100
Clear Non-Vented	092002BBG	100



2.8mm x 19mm	Part No.	Qty/Bag
Solid Vented	091013BBG	100
Solid Non-Vented	091003BBG	100
Tinted Vented	093013BBG	100
Tinted Non-Vented	093003BBG	100
Clear Vented	092013BBG	100
Clear Non-Vented	092003BBG	100



3.2mm x 25.4mm	Part No.	Qty/Bag
Solid Vented	091014BBG	100
Solid Non-Vented	091004BBG	100
Tinted Vented	093014BBG	100
Tinted Non-Vented	093004BBG	100
Clear Vented	092014BBG	100
Clear Non-Vented	092004BBG	100



5mm x 25mm	Part No.	Qty/Bag
Solid Vented	091015BBG	100
Solid Non-Vented	091005BBG	100
Tinted Vented	093015BBG	100
Tinted Non-Vented	093005BBG	100
Clear Vented	092015BBG	100
Clear Non-Vented	092005BBG	100



10mm x 19mm	Part No.	Qty/Bag
Solid Vented	091016BBG	100
Solid Non-Vented	091006BBG	100
Tinted Vented	093016BBG	100
Tinted Non-Vented	093006BBG	100
Clear Vented	092016BBG	100
Clear Non-Vented	092006BBG	100



2mm x 9mm x 25mm	Part No.	Qty/Bag
Solid Vented	091017BBG	100
Solid Non-Vented	091007BBG	100
Tinted Vented	093017BBG	100
Tinted Non-Vented	093007BBG	100
Clear Vented	092017BBG	100
Clear Non-Vented	092007BBG	100



2mm x 16mm x 25mm	Part No.	Qty/Bag
Solid Vented	091018BBG	100
Solid Non-Vented	091008BBG	100
Tinted Vented	093018BBG	100
Tinted Non-Vented	093008BBG	100
Clear Vented	092018BBG	100
Clear Non-Vented	092008BBG	100



3mm x 25mm x 25mm	Part No.	Qty/Bag
Solid Vented	091019BBG	100
Solid Non-Vented	091009BBG	100
Tinted Vented	093019BBG	100
Tinted Non-Vented	093009BBG	100
Clear Vented	092019BBG	100
Clear Non-Vented	092009BBG	100

Assorted – Solid	Part No.	Qty/Bag
Solid Vented	091020BBG	100
Solid Non-Vented	091010BBG	100

Assorted – Tinted	Part No.	Qty/Bag
Tinted Vented	093020BBG	100
Tinted Non-Vented	093010BBG	100

Assorted – Clear	Part No.	Qty/Bag
Clear Vented	092020BBG	100
Clear Non-Vented	092010BBG	100

DuoGuard™ Instrument Tip Protectors

- Ensure that your instruments stay in the proper open position during processing with DuoGuard™ Instrument Tip Protectors
- Available vented and non-vented, they are tinted for easy viewing of instrument tips
- Not made with natural rubber latex and intended for single-use




Tinted

Part No. Vented	Part No. Non-Vented	Color	Size mm	Qty
093031bbg	093021BBG	White	6.4 x 19	50
093032bbg	093022BBG	Blue	9.5 x 20	50
093033bbg	093023BBG	Green	12.7 x 20	50
093034bbg	093024BBG	Red	19 x 25.5	50
093035bbg	093025BBG	Yellow	15 x 25	50
093037bbg	093027BBG	Orange	25 x 25	50
093040bbg		Assorted	Assorted	50

Low Temperature Web Guards

- An Aspen exclusive
- These unique guards with a flow through, open weave are designed to be used only with low heat, plasma sterilization processes, such as the Sterrad™ Sterilization System
- Great for helping to protect endoscopes and their accessories
- Not made with natural rubber latex and intended for single-use




Part No.	Color	mm	Qty/Bag
099001BBG	Orange Mesh	10 x 76	24
099002BBG	Green Mesh	15 x 76	24
099003BBG	Blue Mesh	20 x 76	24
099004BBG	Red Mesh	40 x 76	24
099005BBG	Black Mesh	60 x 76	16

InnerVent™ Instrument Tip Protectors

- Offer the performance of vented tip protectors with the safety of non-vented tip protectors
- Raised inner nubs allow improved flow through of steam and gas while preventing the tip protectors from “vacuuming” onto the instrument during drying
- Not made with natural rubber latex and intended for single-use



Tinted



Solid/Radiopaque



Clear



Solid, Clear, Tinted

Part No.	Color	Size mm	Qty/Bag
091117BBG	Orange, Solid	2 x 9 x 25	100
091118BBG	Brown, Solid	2 x 16 x 25	100
091119BBG	Black, Solid	3 x 25 x 25	100
092117BBG	Clear	2 x 9 x 25	100
092118BBG	Clear	2 x 16 x 25	100
092119BBG	Clear	3 x 25 x 25	100
093117BBG	Orange, Tinted	2 x 9 x 25	100
093118BBG	Brown, Tinted	2 x 16 x 25	100
093119BBG	Black, Tinted	3 x 25 x 25	100

Gelpi Retractor Tip Protectors

- Protect expensive Gelpi retractors with Aspen Surgical's exclusive Gelpi Retractor Tip Protectors
- Available in three sizes to fit the most commonly used Gelpi retractors
- Gas and steam sterilization compatible
- Not made with natural rubber latex and intended for single-use



Part No.	Color	Size mm	Qty/Bag
096109BBG	Clear	38 Length	55
096129BBG	Clear	51 Length	45
096139BBG	Clear	76 Length	40

Weitlaner Instrument Tip Protectors

- Help Protect your expensive Weitlaner instruments and reduce reprocessing costs associated with compromised wrap
- Specifically designed for use with up to a 4" Weitlaner retractor and offer easy one-handed application
- Not made with natural rubber latex and intended for single-use



Vented, Solid or Clear

Part No.	Size	Color	Qty/Bag
096409BBG	Small	Magenta	50
096410BBG	Small	Clear	50
096419BBG	Medium	Green	50
096420BBG	Medium	Clear	50

Laparoscopic MIS Tip Protectors

- Venting of these specialty tip protectors provides extra flow through to Trocar probes and Veress needles during high level disinfection or sterilization
- Not made with natural rubber latex and intended for single-use



Vented and Clear

Part No.	Size mm	Applications	Qty/Bag
097011BBG	19	Use for Veress Needle	50
097012BBG	20 x 5	Use for Trocar Probe	25
097013BBG	43 x 5	Use for Trocar Probe	20
097014BBG	20 x 10	Use for Trocar Probe	25
097015BBG	57 x 10	Use for Trocar Probe	20

Endoscope Guards

- Help safeguard the eye pieces of delicate endoscopic equipment during storage with these soft, conformable guards
- Not made with natural rubber latex and intended for single-use

Blue, Non-Vented

Part No.	Description	Qty/Bag
095016BBG	Small Eye Piece	6



Round Tip Protectors

- Designed to help protect valuable instruments, ranging in sizes from the smallest Veress Needle to the largest endoscope, during handling and storage
- Tinted for easy visual inspection and vented for improved steam, gas, flow
- Not made with natural rubber latex and intended for single-use



Vented and Tinted

Part No.	Part No.	Color	Fits Tip Sizes	Qty/Bag
096001BBG	096021BBG	White	2-3	70
096002BBG	096022BBG	Blue	3-4	70
096003BBG	096023BBG	Green	4-5	70
096004BBG	096024BBG	Red	5-6	70
096005BBG	096025BBG	Yellow	6-7	70
096006BBG	096026BBG	Maroon	7-8	70
096007BBG	096027BBG	Orange	8-9	70
096008BBG	096028BBG	Brown	9-10	70
096009BBG	096029BBG	Black	10-11	70
096013BBG	096033BBG	Green	13-14	50
096015BBG	096035BBG	Yellow	15-16	50

ColorTrack® Instrument Identification Tape

- Helps make organizing and tracking valuable instruments easy
- It is flexible, adheres securely to instruments, and will not chip, flake, peel, or crack
- With 75+ colors, symbols, and character options, the color-coding options are nearly endless
- If a change is needed, simply remove the ColorTrack I.D. Tape and replace with your new color choice
- Not made with natural rubber latex and non-sterile

Sheet Tape

- Each 8½" x 11" sheet has 374 total inches of tape and contains:
44 each of 2 inch strips
44 each of 1½ inch strips
220 each of 1 inch strips

Solid Color
¼" Height / Various widths

Part No.	Description
152001EEA	White
152002EEA	Royal Blue
152003EEA	Kelly Green
152004EEA	Red
152005EEA	Yellow
152006EEA	Purple
152007EEA	Orange
152008EEA	Brown
152009EEA	Black
152010EEA	Gray
152011EEA	Lavender
152012EEA	Beige
152013EEA	Light Blue
152014EEA	Lime Green
152015EEA	Pink
152017EEA	Plum
152018EEA	Mint Green

Dots
¼" Height/Variou widths

Part No.	Description
158040EEA	Dots, Fuchsia
158041EEA	Dots, Medical Blue
158042EEA	Dots, Emerald Green
158043EEA	Dots, Violet

Stars
¼" Height / Various widths

Part No.	Description
158030EEA	Stars, Fuchsia
158031EEA	Stars, Medical Blue
158032EEA	Stars, Emerald Green
158033EEA	Stars, Violet

White Stripe
¼" Height / Various widths

Part No.	Description
154002EEA	White/Royal Blue
154003EEA	White/Kelly Green
154004EEA	White/Red
154005EEA	White/Yellow
154007EEA	White/Orange
154008EEA	White/Brown
154009EEA	White/Black
154011EEA	White/Lavender
154013EEA	White/Light Blue
154014EEA	White/Lime Green
154015EEA	White/Pink
154017EEA	White/Plum

Diagonal Stripes
¼" Height/Variou widths

Part No.	Description
159001EEA	Diagonal, Purple/Black
159002EEA	Diagonal, Orange/Black
159003EEA	Diagonal, Blue/Black
159004EEA	Diagonal, Rubine/Black
159005EEA	Diagonal, Yellow/Black
159006EEA	Diagonal, Green/Black
159007EEA	Diagonal, White/Black

Tri-Color Stripe
¼" Height / Various widths

Part No.	Description
155051EEA	Red/White/Blue
155052EEA	Red/White/Black
155053EEA	Red/Yellow/Black
155054EEA	Red/Black/Green

Black Stripe
¼" Height / Various widths

Part No.	Description
153001EEA	Black/White
153003EEA	Black/Kelly Green
153004EEA	Black/Red
153005EEA	Black/Yellow
153007EEA	Black/Orange
153010EEA	Black/Gray
153011EEA	Black/Lavender
153012EEA	Black/Beige
153013EEA	Black/Light Blue
153014EEA	Black/Lime Green
153015EEA	Black/Pink
153017EEA	Black/Plum

Triangles
¼" Height/Variou widths

Part No.	Description
158050EEA	Triangles, Fuchsia
158051EEA	Triangles, Medical Blue
158052EEA	Triangles, Emerald Green
158053EEA	Triangles, Violet

Numbers
¼" Height / Various widths

Part No.	Description
158020EEA	Numbers, Fuchsia
158021EEA	Numbers, Medical Blue
158022EEA	Numbers, Emerald Green
158023EEA	Numbers, Violet

Assorted Symbols
¼" Height / Various widths

Part No.	Description
158060EEA	Red Hearts on White
158080EEA	Alphabet, Fuschia/Wt, A-E

4-Color
¼" Height / Various widths

Note: Assorted color package will contain separate sheets each with its own color.

Part No.	Description
152031EEA	Red/Yellow/Purple/Orange
152032EEA	Blue/White/Black/Green

Specialty Stripe
¼" Height / Various widths

Part No.	Description
156003EEA	Blue/Green
156004EEA	Blue/Red
157002EEA	Yellow/Blue
157003EEA	Yellow/Green
157004EEA	Yellow/Red

Code to Specialty
¼" Height / Various widths

Part No.	Description
ENT	Mint/White ENT
ER ER	Purple/White ER
GEN	Royal Blue/White GEN
OB/GYN	Lt Blue/White OB/GYN
ORTHO	Kelly Green/White ORTHO
NEURO	Orange/White NEURO
CARDIAC	Red/ White CARDIAC
6	Red/White ASPEN SYMBOL

Roll Tape
Roll Tape ¼" x 300"

Part No.	Description
151001EEA	White
151002EEA	Royal Blue
151003EEA	Kelly Green
151004EEA	Red
151005EEA	Yellow
151006EEA	Purple
151007EEA	Orange
151008EEA	Brown
151009EEA	Black

Roll Tape ¼" x 250"

Part No.	Description
151011EEA	White
151012EEA	Royal Blue
151013EEA	Kelly Green
151014EEA	Red
151015EEA	Yellow
151016EEA	Purple
151017EEA	Orange
151018EEA	Brown
151019EEA	Black



Application Instructions:
For secure placement through multiple processing cycles, ColorTrack® I.D. Tape should be wrapped 1½ times around the Instrument and excess length should be trimmed. Instrument sterility under ColorTrack® I.D. Tape is assured after processing with properly applied tape.*

* Test Results Available On Request
Size and color of tapes shown may vary from actual products.

Sterion® Same Day Tape Remover

- An easy, instrument-friendly solution to remove ColorTrack® Tape
- Helps eliminate the risk of cuts and sticks common with scalpel or utility knife tape removal
- Neutral scent helps make Same Day Tape Remover comfortable to use


Part No.	Description
151066EEA	1 Gallon Container

Helps to
Easily
Remove
Our Tape!



ColorTrack® Write-On ID Labels


- Ideal for customization of trays, storage areas, cabinets or carts, these color-coded labels are formatted on 8½”x 11” sheets with three different label sizes on each sheet
- For best results use with permanent markers only
- Not made with natural rubber latex



Part No.	Description	Qty/Sheet
161002EEA	Blue Border	44
161003EEA	Green Border	44
161004EEA	Red Border	44
161005EEA	Yellow Border	44

Oxboro® Corner Protectors

- Protect your instrument trays from rips, holes, and tears in sterilization wrap due to sharp tray corners, improper handling, or heavy surgical trays
- Aspen Surgical’s Corner Protectors are designed to create a flexible barrier between the tray’s sharp corners and the sterilization wrap
- They can be used with steam, gas, ozone, and low temperature plasma sterilization and are reusable
- Not made with natural rubber latex




Part No.	Description	Qty/Bag
CP-80NS	Corner Protectors	80

Sterion® Color-Coded Identification Locking Tags

- These tamper evident locking tags can be customized to a variety of uses such as identifying instruments, carts, trays, containers or cabinets
- Available pre-printed with “Repair”, “Sharpen”, or “Good” on one side and writable on the other or blank on one side and writable on the other
- Not made with natural rubber latex



Part No.	Description	Qty/Bag
251001BBG	Red, Imprint: Repair, Write-On	100
251002BBG	Yellow, Imprint: Sharpen, Write-On	100
251003BBG	Green, Imprint: Good, Write-On	100
251004BBG	Blue, Write-On	100
251005BBG	White, Write-On	100
251006BBG	Green, Write-On	100
251007BBG	Black, Write-On	100
251008BBG	Purple, Write-On	100
251009BBG	Orange, Write-On	100
251010BBG	Assorted Colors, Write-On	100




SafeClean® Instrument and Channel Cleaning Brushes

- Designed to clean surgical/medical instruments and endoscopic devices
- The complete line of SafeClean® Cleaning Brushes are strong enough to effectively remove surgical residue, yet gentle enough to use on the most delicate instruments and devices
- Not made with natural rubber latex

General Purpose Brushes

- Can be used throughout all of your instrument and device reprocessing areas to clean uniquely shaped and larger bore surgical, medical, and endoscopic instruments and devices
- Not made with natural rubber latex



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Qty/ Bag
241007BBG	1mm (0.039")	18"	0.5"	2
241010BBG	1.6mm (0.063")	18"	0.5"	2
241015BBG	1.8mm (0.071")	18"	0.5"	2
241019BBG	2mm (0.079")	12"	1"	2
241020BBG	2.2mm (0.086")	18"	0.5"	2
241024BBG	2.5mm (0.098")	24"	2.5"	2
241025BBG	2.7mm (0.106")	18"	1"	2
241029BBG	3mm (0.118")	12"	1"	2
241030BBG	3.2mm (0.126")	18"	2"	2
241032BBG	4mm (0.158")	12"	1"	2
241033BBG	5mm (0.197")	12"	2.5"	2
241034BBG	5mm (0.197")	24"	2.5"	2
241035BBG	5mm (0.197")	16"	2.5"	2
241040BBG	6mm (0.236")	16"	2.5"	2
241041BBG	6mm (0.236")	24"	2.5"	2
241042BBG	7mm (0.275")	16"	2.5"	2
241043BBG	7mm (0.275")	24"	2.5"	2
241044BBG	10mm (0.393")	16"	2.5"	2
241045BBG	11mm (0.433")	16"	2.5"	2
241046BBG	11mm (0.433")	24"	2.5"	2
241047BBG	10mm (0.393")	24"	2.5"	2
241050BBG	12mm (0.472")	16"	2.5"	2
241055BBG	15mm (0.590")	16"	2.5"	2
241060BBG	20mm (0.787")	16"	2.5"	2

Kits

- Help make ordering easier with all of the brushes you need in one kit

Small Cleaning Brush Kit

- Includes various small general purpose cleaning brushes
- Not made with natural rubber latex

Part No.
241100bbg – 1 Each
241010BBG
241015BBG
241020BBG
241025BBG
241030BBG

Large Cleaning Brush Kit

- Includes various large general purpose cleaning brushes
- Not made with natural rubber latex

Part No.
241200bbg – 1 Each
241035BBG
241040BBG
241045BBG
241050BBG
241055BBG
241060BBG

“Best Sellers” Cleaning Brush Kit

- Includes a selection of our most popular general purpose cleaning brushes
- Not made with natural rubber latex

Part No.
241300bbg – 1 Each
241010BBG
241015BBG
241020BBG
241025BBG
241030BBG
241035BBG
241040BBG
241045BBG
241050BBG
241055BBG
241060BBG

General Cleaning Brush Kit

- Contains an extensive selection of general purpose brushes, as well as toothbrush and fan tip brushes
- Not made with natural rubber latex

Part No.
241600bbg – 3 Each
241025BBG
241030BBG
241035BBG
241040BBG
241041BBG
241045BBG
241046BBG
241050BBG
241055BBG
241060BBG
242001BBG
242002BBG
242003BBG
243001BBG
243004BBG

Micro/Ophthalmic Brush Kit

- Contains a selection of our smallest general purpose cleaning brushes
- Not made with natural rubber latex

Part No.
241700bbg – 3 Each
241007BBG
241010BBG
241015BBG
241020BBG
241025BBG
241030BBG


Tooth Brush Kit

- Not made with natural rubber latex

Part No.
243100bbg – 1 Each
243001BBG
243002BBG
243003BBG

Suction Tube Cleaning Brushes

- Brushes are designed to fit into and effectively clean very small bore devices and instruments used in surgery, MIS, and therapy procedures today
- Not made with natural rubber latex



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Qty/ Bag
241003BBG	3 French (1mm, 0.039")	6"	0.5"	2
241005BBG	5 French (1.67mm, 0.066")	6"	0.5"	2
241006BBG	6 French (2mm, 0.079")	8"	0.5"	2
241008BBG	7 French (2.3mm, 0.092")	8"	0.5"	2
241009BBG	8 French (2.7mm, 0.105")	8"	0.5"	2
241011BBG	9 French (3mm, 0.118")	8"	0.5"	2
241012BBG	10 French (3.3mm, 0.131")	8"	0.5"	2
241013BBG	11 French (3.7mm, 0.144")	8"	0.5"	2
241014BBG	12 French (4mm, 0.158")	8"	0.5"	2

Fan-Tip Brush

- Features bristles that wrap around the end of the brush to clean bottles, jars, and other large instruments
- Not made with natural rubber latex



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Qty/ Bag
242001BBG	4.8mm-12.7mm Tapered (0.189"–0.5")	7"	3"	2
242002BBG	13mm (0.512")	6"	3"	2
242003BBG	16mm (0.623")	5"	2.5"	2
242046BBG	9.6mm (0.384")	24"	2.5"	2

Tracheostomy Brushes


- Uniquely designed to clean tubes
- Not made with natural rubber latex



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Qty/ Bag
247001BBG	Fits Trach Sizes 00–5FR 9mm (0.354")	4.75"	1.85"	2
247002BBG	Fits Trach Sizes 6–12FR 13mm (0.512")	8"	3.66"	2

Toothbrush Style

- Variety of bristle types available
- Designed to clean all types of instruments and stains including nylon, brass, and steel
- Not made with natural rubber latex



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Qty/ Bag
243001BBG	13mm (0.512"), Nylon	7.25"	1.5"	2
243002BBG	13mm (0.512"), Brass	7.25"	1.5"	2
243003BBG	13mm (0.512"), Steel	7.25"	1.5"	2

Double-Ended Nylon Cleaning Brushes

- Features one end with an extra large brush head for general cleaning and one end with a single row of bristles for more intricate cleaning
- Not made with natural rubber latex




Available Sterile and Non-Sterile

Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Qty
243004BBG	Non-Sterile			
	Short End, 2mm (0.079")	7"	0.867"	2/Bag
	Long End, 8mm (0.315")		1.38"	
243014PBX	Sterile			
	Short End, 2mm (0.079")	7"	0.867"	1/Pack,
	Long End, 8mm (0.315")		1.38"	12/Packs/Box

Reusable Channel Cleaning Brushes


- Channel Brushes with tapered brush heads, elevator sheath, and comfortable round, plastic handles
- Available in a variety of lengths and materials to fit the range of endoscopic devices
- Not made with natural rubber latex



Part No.	Description	Total Length	Sheath Diameter	Qty/ Bag
248107BBG	3-5mm	70cm	1.6mm	1
248115BBG	3-5mm	115cm	1.6mm	1
248160BBG	3-5mm	160cm	1.6mm	1
248180BBG	3-5mm	230cm	1.8mm	1
248190BBG	3-5mm	120cm	1.8mm	1

Disposable Scope Channel Cleaning Brushes


- Brush heads that range from to 2.5mm - 5mm and will fit scopes as small as 1.1mm
- Features more than double the bristle length than competing low-cost brushes
- Strong, sturdy handle making it easy to pull through the scope, while still cleaning the scope effectively
- Intended for single-use and not made with natural rubber latex



Part No.	Description	Total Length	Sheath Diameter	Qty/ Box
248102BBG	2.5mm	230cm	1mm	50
248103BBG	3mm	230cm	1.55mm	50
248104BBG	5mm – Double Ended	230cm	1.85mm	50
248105BBG	5mm	230cm	1.85mm	50

Hand/Nail Brushes

- Used to clean dirt on hands and nails
- Can also be used for general instrument cleaning
- Not made with natural rubber latex



Part No.	Brush Head Area	Qty/Bag
244001BBG	35mm x 93mm (1.375" x 1.75")	2

Instrument Brush Storage

- Convenient storage with our acrylic racks organize your cleaning brushes
- Helps ensure you will have the right brush for the job
- Not made with natural rubber latex

Part No.	Description	Qty/Box
246001BBG	Brush Rack, Medium (without drip tray)	1
246003BBG	Brush Rack, Large (with drip tray)	1



Operating Room Disposables

Disposable Camera Covers and Light Handles

Helps save time and reduce inconvenience associated with having to sterilize reusable light handles and camera covers. Adaptors are currently available for many common O.R. lights and cameras.

O.R. Camera Cover

Part No.	Description	Qty/Box
LT-C02	Camera Cover	20



Specialty Disposable Light Handle

Part No.	Description	Sterility	Qty/Box
LT-ALC01	Disposable Light Handle	Sterile	18




Standard Disposable Light Handle

Part No.	Description	Qty/Box
LT-R01	Rigid Handle	24
LT-F01B	Flexible Cover (1 Each)	100
LT-F02	Flexible Cover (2 Each)	80
LT-F03	Flexible Cover (3 Each)	80



SurgiSafe® Absorbent Floor Mats

Available in multiple sizes and absorbency levels, the absorbent floor mats help keep floors clean, dry, and slip-free, protecting your staff from slips, trips, and falls.



Part No.	Color	Absorbency Level	Fluid Barrier Backing	Size	Qty/Case
82350	Blue	Medium	Yes	23" x 40"	50
84625	Blue	Medium	Yes	46" x 40"	25
84610	Blue	Medium	Yes	46" x 40"	10
84672	Blue	Medium	Yes	40" x 72"	10
81860-2	Yellow	Medium	No	18" x 40"	60
83630-2	Yellow	Medium	No	36" x 40"	30
83610-2	Yellow	Medium	No	36" x 40"	10
83615-2L	Yellow	Medium	No	36" x 72"	15
83625-2L	Yellow	Medium	No	36" x 72"	12
83630-M	Yellow	Medium	No	36" x 40"	15
	White	Medium	No	28" x 40"	15

ErgoSupport® Anti-Fatigue Mat

Available in 5 sizes and water-resistant, this anti-fatigue mat is intended to help reduce fatigue that healthcare personnel experience during long hours of standing in one location.



Part No.	Description	Color	Size	Qty/Case
81005	Anti-Fatigue Mat	Red	13" x 15"	5
82001	Anti-Fatigue Mat	Red	18" x 24"	1
82005	Anti-Fatigue Mat	Red	18" x 24"	5
83002	Anti-Fatigue Mat	Red	20" x 48"	2
84003	Anti-Fatigue Mat	Red	20" x 32"	3

Hill-Rom is a leading global medical technology company with more than 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Around the world, Hill-Rom's people, products, and programs work towards one mission: **Every day, around the world, we enhance outcomes for patients and their caregivers.**

Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2016 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.

M-SM-CAT-1020821 rev 1 18-APR-2016 ENG – US

For further information about this product or a service, please contact your local Aspen Surgical representative or visit our webpage:

USA 616-698-7100

Canada 888-364-7004

www.aspensurgical.com

Enhancing outcomes for
patients and their caregivers:

Hill-Rom