



The *smarter source* for affordable specialty instruments

## MASTEL TWO-RAY CORNEAL MERIDIAN MARKER

### INTRA-OP MARKERS

Bores Handle



Thornton Handle



- Designed for Toric IOL's, LRI and AK markings
- 9.0mm Inside Diameter to 11.5mm Outside Diameter Fins
- Fins designed to follow the curvature of the eye
- Fins are sharpened to impress/divide the epithelium without tearing


**Bores**  
**MRMB-2R**

**Thornton**  
**MRMT-2R**



Ambler Surgical is an authorized distributor for the Mastel Precision line of ophthalmic instruments and diamond scalpels. Mastel Precision has been involved exclusively in the design and production of ophthalmic instruments since 1981. Mastel Precision offers unparalleled design, quality and precision-crafted surgical instruments and diamond scalpels for astigmatism management and corneal procedures that are considered unmatched in the ophthalmic industry.

To view pricing and our complete line of Mastel Precision Instruments, visit us online.



OPHTHALMOLOGY

[amblersurgical.com/mastel](http://amblersurgical.com/mastel)

1-888-407-0006