

## Sterilization Statement

### ID Tape

The Aspen Surgical Identification Tapes are intended to identify surgical instruments through processing. They have been evaluated for the following sterilization methods:

Conventional Steam Autoclave Cycle for 121 degrees Celsius for 30 minutes


Steam Autoclave Cycle at 132 degrees Celsius for at least 3 minutes

Conventional Hospital Ethylene Oxide Cycle

After each of the above cycles, all test samples demonstrated to be sterile. The tapes are gas and steam permeable.

Caledonia, Michigan 02DEC2022

Place and Date of Issue

A handwritten signature in black ink that reads "Kathleen Herremans".

Kathleen Herremans  
Kathleen Herremans (Dec 2, 2022 05:21 PST)  
Regulatory Compliance Specialist  
Aspen Surgical Products, Inc.