



## MEDIN STAINLESS STEEL TRAYS

### MEDIN Stainless Steel Trays

MEDIN instrument trays are designed to encase instruments during processing and are safe for all standard sterilization methods: Autoclave/dry heat up to 360°F, ETO, and Cold Solutions.

- An appropriate sterile barrier will be required to maintain sterility.
- Use of the Medin Instrument Trays must be in accordance with the operating instruction supplied with your particular sterilization equipment and in keeping with your particular sterilization equipment and in keeping with your hospital policy on sterilization validation. Several organizations (e.g., AAMI, AORN and ISO) publish guidance information on sterilization using should need additional information.
- Periodic cleaning of the trays and silicone mats/bars after sterilization process is recommended to prolong usable life.
- Do not use if damaged.